



MVP – VACUUM STOPPERING SYSTEM

The Cozzoli MVP provides reliable, semi-automatic, batch insertion of stoppers, plungers or plugs into filled containers. This high performance machine is the perfect companion to the Cozzoli F329 Series liquid filling machines.

A pneumatic stopper transfer mechanism is operated by the cycle control unit. This cycle control unit governs all functions of the assembly operation according to operator selected parameters. Unlike other systems, the vacuum chamber and all container handling parts are completely autoclavable, designed to withstand temperatures of 450° F (232° C). The entire 316L stainless steel chamber can be steam sterilized for maximum cleanliness.

◀ www.cozzoli.com ▶

50 Schoolhouse Rd., Somerset, NJ 08873 | 732.564.0400 | Sales@Cozzoli.com



OPERATING PRINCIPLE

The operator places the filled syringes into a rigid stainless steel barrel rack, tops it with a transfer plate loaded with pre-packaged or individual stoppers, then places the assembly into the MVP's vacuum chamber. The MVP's control unit then regulates the user's vacuum, nitrogen and compresses air utilities which are programmed to remove entrapped air from the product, purge the chamber and transfer the stoppers into the syringes. Final venting to ambient air pressure forces the stoppers into contact with the liquid for a secure, reliable seal.

CONTROL FEATURES

Vacuum Only and Vacuum/Purge/Vacuum Switches:

By setting selector switches on the front panel, the operator can choose to run the system using vacuum only, or with a vacuum/nitrogen purge/vacuum sequence.

Vacuum Level Settings: Internal adjustments allow the user to preset the necessary level of vacuum. The rest of the stoppering cycle will not activate until this level is attained.

Purge Interval (Preset Internally): When Vacuum/Purge/Vacuum is selected and the initial vacuum state is complete, nitrogen automatically feeds into the chamber for the pre-set interval to displace any possible remaining gasses. A low-nitrogen detection feature automatically stops the cycle and warns if nitrogen pressure is inadequate.

Second Vacuum Interval Timer: A digital timer on the front panel sets the interval necessary for complete stopper transfer.

Stoppering Speed: This internal adjustment allows independent control of the advance/retract action of the transfer punch-plate.

Stoppering Force: The force of the pneumatic transfer punch-plate that transfers stoppers into the containers is monitored and controlled by an air pressure gauge and regulator on the front panel.

Venting: Once set, all of the above vacuum, purge and stopper transfer function are automatic. When completed, the system automatically vents to ambient air pressure and the stoppers are forced into final seating position. The rack of containers is then ready for removal.

STANDARD FEATURES


- 316 Stainless Steel Vacuum Chamber
- 115V, 60Hz, 1 Ph Electricals
- Transfer Punch-Plate Assembly
- Available Patterns 5x6, 6x7, 8x8, 10x10 or 10x16
- Stopper Transfer Plate
- Available Patters 5x6, 6x7, 8x8, 10x10 or 10x16
- Container Handling Rack
- Available Patterns 5x6, 6x7, 8x8, 10x10 or 10x16

OPTIONAL FEATURES

- Additional Transfer Punch-Plate Assembly
- Additional Transfer Stopper Transfer Plates
- Additional Container Handling Rack
- Electrical Configuration other than 115V, 60HZ, 1Ph
- MVP-ST with Smaller, Taller Vacuum Chamber
- Dimensions: 6.5"W x 4.75"D x 14.25H ID
- MVP-DT with Narrower, Deeper Vacuum Chamber
- Dimensions: 6.5"W x 17.25"D x 14.25"H ID
- Dual Stage Vacuum Pump
- Exhaust Filter for Vacuum Pump
- Replacement Element for Exhaust Filter
- MVP-CNTL – Controller only
- MVP-STCH – ST Vacuum Chamber Only, No Controller
- MVP-DTCH – DT Vacuum Chamber Only, No Controller
- Spare Parts Kit
- IQ/OQ Validation Protocol
- Factory Acceptance Test Protocol
- PQ Performance Qualification Protocol

REQUIRED SERVICES

- **Electrical:** Grounded 3-wire
- **Standard:** 110-120V, 60 Hz, 5A, 1 Ph
- **Optional:** 220-240V, 50 Hz, 3A, 1 Ph
100-110V, 50 Hz, 6A. 1 Ph



Pneumatic Supply: 60 PSIG (4.2Kg/cm²), .25" (6.4mm) Inside Diameter supply line connects to control unit.

Air Consumption: .04 cu. Ft. (1.2 liter) per cycle at 60 PSIG (4.2kg/cm²) connects to control unit.

Vacuum Supply: In-plant system providing vacuum at 28" (71 cm) Hg. (If in-plant vacuum is not available, a 500 liter/minute vacuum pump will suffice for most applications.)

Nitrogen Supply: Regulated supply for the nitrogen purge feature. (Maximum pressure limited to 3 PSIG (0.2 kg/cm²))

Sterile Filters: Provided by the user.