



## Case History

**Customer:** MedPharmex

**Location:** Pomona, California

**Machine:** GW24

**Application:** Animal Health - MedPharmex is washing large plastic and glass vials that will be subsequently filled with injectable veterinary products. This Cozzoli GW24 is washing vials up to 1000 ml.



### Features

Customized Washing Cycles

Low Maintenance

Easy Changeover

Sanitary Valves and  
Supply Connections

Handles Various Wash Media  
Including Sterile Air, Water &  
Steam

The GW 24 Batch Washer operation simply requires placing a loaded box plate of containers into the washer, positioning it on the center plate, lowering the cover, and push start. Air, steam and 3 different liquids may be used for internal washing, while adjustable spray heads wash container exteriors.



## Case History, Continued



### *More Features*

Accurate,  
Repeatable Washing Cycles

Modified with a Higher Lid for  
Tall Bottles

Provided with a Complete IQ/OQ  
Documentation Package

Allen Bradley SLC 500 Series PLC  
Provides for Storage of up to  
9 Wash Recipes of 18 Steps Each

MedPharmex expanded their production with this GW24. A competitor of theirs was shut down by the Food & Drug Administration and MedPharmex had an opportunity to add to their customer base. Choosing the Cozzoli Washer confirmed their commitment to quality with their new customers, and allowed them to quickly take advantage of this unexpected opportunity due to the typically fast delivery on this model.

*For more information on this filler, contact the Cozzoli Machine Company Sales Department phone: 732-564-0400, fax 732-564-0444, email [sales@cozzoli.com](mailto:sales@cozzoli.com), or visit our website: [www.cozzoli.com](http://www.cozzoli.com).*